## CYCLOPHOTOCOAGULATION PROBE STERILE, SINGLE USE

## ENGLISH

### INDICATIONS:

U.S. - This product, when used with an infrared laser, is indicated for Transscleral Cyclophotocoagulation (TSCPC) of the ciliary processes.

Outside U.S. - This product, when used with an infrared laser, is indicated for Transscleral Cyclophotocoagulation (TSCPC) of the ciliary processes for treatment of glaucoma

### DIRECTIONS FOR USE:

#### Connection of probe

A probe with catalog number 620231 is compatible with LIGHTLas 810 and TruScan Pro 810.

1. Release the probe from pouches and packaging and remove protective cap.
2. Connect the probe to a compatible infrared laser machine. Follow the recommendations of the laser source manufacturers (e.g.: calibration of the laser source as well as the need for protective eye wear for the corresponding wavelength of the laser model used).
3. If using the probe in LIGHTMED SP-Mode™, unscrew the footplate carefully from probe exposing only the metal threaded rod and ball fiber.
4. The probe is now ready to use.

If not handled properly during any step of setup (such as bending, improperly securing, poorly inserting, etc.), the glass optical fiber within the product can be damaged leading to harm to the patient or laser operator.

#### Disposal of probe

Follow standard facility procedures for handling of biohazardous material after each use of the delivery device.

#### Anesthesia

Local anesthetic is necessary: retrobulbar or peribulbar injection, or subconjunctival anesthesia, with, for example, 2% mepivacaine, or equivalent agent

#### Eye Moistness

Keep the tip of the probe moist and in contact with the eye, by applying a film of methylcellulose solution to the fiber-optic tip, or ask the patient to close their eyelids, pretreatment, to assure that a natural tear film exists. If you are planning on using a lid speculum to hold the lids open, then please apply artificial tears to the eye, with the speculum in place. Whatever method that you choose to use, please reapply after every four laser applications.

#### Laser Settings

Laser setting selection and adjustment are the responsibility of the clinician, for each patient. General Power settings are 1250 mW for 4000 milliseconds for a patient with a Dark Brown Iris. 1500 mW for 3500 milliseconds for all other iris colors. These settings are effective for most eyes. If unsure what settings are required, select a low power. Using of the parameters laid out in this section will typically result in few or no audible “pops.” Inspect the tip of the optic fiber after each application to verify that it is not visibly blackened, it is free of adherent tissue, it retains its spherical shape, and is not cracked. The aiming beam should emerge from the optic fiber tip with a diverging beam profile and produce a round spot on a perpendicular target with uniform brightness.

### WARNING:

Excess treatment power may result in ocular surface burns or ciliary body hemorrhage. The maximum possible laser power for this device is the lesser of the maximum power of the laser system and 2500 mW. Contamination of the tip or reuse of the probe may result in these types of burns. Excessive energy settings on the laser may result in burn at the ocular equator. Heavy pigmentation at the perilimbal conjunctiva may also result in burns, therefore these areas should be avoided.

This device should be used with an appropriate laser-safety filter. No one should ever look directly into the laser light source or the tip of the probe. All OR personnel should be instructed to wear safety goggles and avoid looking at surfaces that the laser light may be scattered from.

### Instructions for use of probe with footplate attached

#### Laser Settings

Laser setting selection and adjustment are the responsibility of the clinician, for each patient. General Power settings are 1250 mW for 4000 milliseconds for a patient with a Dark Brown Iris. 1500 mW for 3500 milliseconds for all other iris colors. These settings are effective for most eyes. If unsure what settings are required, select a low power. Using of the parameters laid out in this section will typically result in few or no audible “pops.” Inspect the tip of the optic fiber after each application to verify that it is not visibly blackened, it is free of adherent tissue, it retains its spherical shape, and is not cracked. The aiming beam should emerge from the optic fiber tip with a diverging beam profile and produce a round spot on a perpendicular target with uniform brightness.

#### Placement

See “Caution” section if operating on a phakic eye

Figure 1

Hold the probe parallel to the visual axis of the eye and perpendicular to the iris plane. The shorter edge of the footplate of the probe should be held firmly against the area located between the middle and the anterior border of the limbus. Laser delivery is transscleral.

#### Applications

Use one-half of the footplate width to determine the correct spacing between treatment sites. Align the side of the probe over the indented center of the previous treatment.

Figure 2

#### Treatment

This should result in 6-7 laser applications per quadrant, 18-21 applications total, avoiding application in the temporal quadrant

Figure 3

#### Patient Response:

Many surgeons prescribe topical cycloplegics and corticosteroids in anticipation of secondary post-operative inflammation and discomfort.

#### Patient Retreatment, (if necessary):

Begin treatment 45 degrees from the initial treatment. The second 270 degree treatment will cover a half of the untreated quadrant, plus two and a half quadrants from the previously treated area

### Instructions for use of probe WITHOUT footplate:

Figure 4

#### Mark the Eye

Using an ophthalmic marker and an ophthalmic measurement tool, mark the eye approximately 3mm from the limbus on both the superior quadrant and inferior quadrant of the eye like shown below:

Figure 5

#### Placement

Hold the probe approximately perpendicular to the surface of the eye, see Figure 5.

#### Applications and Treatment

Figure 6

To treat the superior quadrant, apply firm pressure and move the device in a continuous arc motion from the 9:30 to 2:30 clock positions. Apply laser continuously, “PAINT”, for 4-7 complete passes between clock positions, for a total of 50-90 seconds in SP-Mode™. Based on recommendations of experienced clinicians: use on the TruScan Pro 810 or LIGHTLas 810 machine, SP-Mode™, at 31.3% (0.5ms on 1.1ms off). To treat the inferior quadrant, repeat process for the suggested time over 3:30 to 8:30 clock hours.

### CYCLOPHOTOCOAGULATION PROBE CLEANLINESS:

Keep the probe tip clean to lower the risk of burns to the ocular surface. If the tip becomes dirty during the procedure, clean it gently with an alcohol swab. Solvents of flammable solutions used for cleaning should be allowed to evaporate before laser equipment is used. If discoloration on the tip cannot be removed by gentle cleaning, discard the Cyclophotocoagulation Probe. Scleral burns are not typical when using this probe. If a burn occurs, discard the probe (it is probably contaminated), and replace it with a new probe. This product is a single use product.

### CAUTION:

This device is to be used only by qualified ophthalmologists after considering patient conditions. After procedure, monitoring the IOP may help prevent secondary post-operative complications. For phakic eyes, the probe must be placed at more than 1mm behind the limbus as to avoid the risk of damage to the lens with the laser beam. In order for correct treatment, the probe must be directly parallel to the visual axis for each laser treatment.

Keep the probe and eye surface moist throughout the treatment.

Do not treat over the site of previous trabeculectomies. If using in SP-Mode™, reduce the clock hours to avoid the previously mentioned while still delivering a total of 50s over the reduced clock hour range.

**PRECAUTION:**

Do not use if package is open, damaged or wet.

Do not use if expiration date shown on the label is exceeded.

**STORAGE:**

Store in a dry area

**ADVERSE REACTIONS:**

Adverse reactions associated with the procedure being performed (to be assessed by the physician) are rare but can include: Ocular surface burns, damage to the lens, phthisis, uveitis, neurotrophic corneal defects, choroidal decompensation, malignant glaucoma, choroidal detachment, atonic pupil, corneal edema, iritis, ophthalmia, necrotizing scleritis, chronic pain, ciliary body hemorrhage, hypotony, IOP spikes, and a decrease in or loss of vision.

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**Trademarks are property of their owners**
Katalyst Surgical, LLC 754 Goddard Ave. Chesterfield, MO USA 63005