



# YAG/SLT/Vitreolysis Quick Start Guide

## YAG Posterior Capsulotomy

YAG Posterior Capsulotomy is the most successful and frequent treatment for secondary cataract.

### Treatment Guidelines

1. Anesthetize the dilated eye with a topical anesthetic
2. A YAG capsulotomy laser lens may help maximize laser focus and stabilize the eye
3. Focus the aiming beam on the posterior capsule in an area that is off the central axis while titrating energy level (going too peripheral may destabilize the IOL)
4. Utilize the lowest possible effective energy level to minimize side effects
5. Recommend no more than 50 mJ per laser session

### Post Treatment

- Fragments of the capsule that are cleared from the central axis are commonly dispersed into the vitreous and may cause temporary blurred vision that typically improves after 24 hours
- Subsequent laser treatments may be considered if the central axis is not adequately cleared
- Patient IOP should be measured 30 minutes after procedure

### Suggested Setting on LIGHTLas YAG: Capsulotomy

|                             |   |
|-----------------------------|---|
| Laser Offset                | +150 posterior offset (to avoid lens pitting) |
| Energy (with laser lens)    | 1.0 - 2.5 mJ                                  |
| Energy (without laser lens) | 0.9 - 2.0 mJ                                  |
| Total Treatment             | Max of 50 mJ per laser session                |

## YAG Peripheral Iridotomy (PI)

YAG Iridotomy treats Narrow Angle Glaucoma by creating an opening in the peripheral iris to release outflow of aqueous humor and reduce the IOP.

### Treatment Guidelines

1. Place a YAG Iridotomy laser lens on the surface of the anesthetized cornea in an undilated patient
2. Focus the YAG laser beam on the peripheral iris to create a small opening
3. Recommend no more than 30 mJ per laser session

### Post Treatment

- Patient IOP should be measured 30 minutes after procedure

### Suggested Setting on LIGHTLas YAG: Iridotomy

|                     |                                |
|---------------------|--------------------------------|
| Laser Offset        | 0 or +150 posterior offset     |
| Energy (light iris) | 2.5 - 3.5 mJ                   |
| Energy (dark iris)  | 3.5 - 4.5 mJ                   |
| Total Treatment     | Max of 30 mJ per laser session |



1130 Calle Cordillera | San Clemente, CA 92673 | USA  
T: 949-218-9555 | F: 949-218-9556 | sales@lightmed.com  
[www.lightmed.com](http://www.lightmed.com)

# YAG/SLT/Vitreolysis

## Quick Start Guide

### YAG Vitreolysis

YAG laser vitreolysis is a highly effective outpatient procedure that can eliminate floaters and vitreous opacities.

#### Treatment Guidelines

1. Recommend a mid-vitreous lens to focus and clearly visualize the floater
2. Ensure the floater is in focus, but the retina must be out of focus to verify adequate distance for treatment (at least 3mm from retina)
3. Commence with low energy and gradually increase the power until there is adequate vaporization
4. Recommend no more than 1000 mJ per laser session

#### Post Treatment

- Patients may see dark specks in their lower field of vision in the initial 15-30 minutes post-procedure, corresponding to micro gas bubbles, which should quickly dissipate
- Patient IOP should be measured 30 minutes after procedure

#### Suggested Setting on LIGHTLas YAG-V: Vitreolysis

|                          |                                  |
|--------------------------|----------------------------------|
| Laser Offset             | 0                                |
| Energy (with laser lens) | 2.5 – 4.5 mJ                     |
| Total Treatment          | Max of 1000 mJ per laser session |

### Selective Laser Trabeculoplasty (SLT)

SLT provides a highly effective and minimally invasive treatment for glaucoma patients in reducing intraocular pressure (IOP).

#### Treatment Guidelines

1. Recommend utilization of an SLT gonio laser lens without magnification to avoid changes to the spot size
2. Laser should be focused on the trabecular meshwork (TM)
3. Titrate the energy at 0.1 mJ increments until the appearance of micro-cavitation bubbles are observed next to the TM
4. Treat area with contiguous, non-overlapping spots along the TM

#### Post Treatment

- Follow up visits are typically scheduled at 1 week, 4 weeks, and 3 months after treatment
- Patient IOP should be measured 30 minutes after procedure

#### Suggested Setting on LIGHTLas SLT

|   |  |
|---|--|
| Energy (lighter TM pigment may require higher energy) | 0.4 - 0.8 mJ                                     |
| Spot Size   | 400 $\mu$ m (standard setting)                   |
| Pulse Width   | 3 ns (standard setting)                          |
| Total Treatment                                       | ~ 50 shots per 180° / or<br>~ 100 shots per 360° |



1130 Calle Cordillera | San Clemente, CA 92673 | USA  
T: 949-218-9555 | F: 949-218-9556 | sales@lightmed.com  
[www.lightmed.com](http://www.lightmed.com)