



LIGHTMED is proud to distribute and represent high-quality laser probes that are FDA and CE approved. This latest design of fiberoptic delivery systems include new flexible curved designs and illuminated probes manufactured with a new design to maximize illumination. These probes are designed to connect to the following lasers: LIGHTLas 532, LIGHTLas 577, LIGHTLas 810, TruScan, and TruScan Pro models.

## Indications for Use:

### Fiberoptic Probes

**20, 23, 25 or 27 gauge, sterile, single use**  
(0.81mm, 0.64mm, 0.50mm, 0.41mm)

#### Description/Indication for Use *(by physicians in the operating room):*

Laser Probes are intended for use in vitreoretinal surgery to perform endo-ocular laser photocoagulation treatments at operating wavelengths of 500nm to 900nm.

Illuminated Laser Probes are intended for use in vitreoretinal surgery to perform endo-ocular laser photocoagulation treatments at operating wavelengths of 500nm to 900nm and to illuminate a surgical site.

### Cyclophotocoagulation Probe

**Sterile, single use**

#### Indications:

**U.S.** - This product, when used with an infrared laser, is indicated for Transscleral Cyclophotocoagulation (TSCPC) of the ciliary processes.

**Outside U.S.** - This product, when used with an infrared laser, is indicated for Transscleral Cyclophotocoagulation (TSCPC) of the ciliary processes for treatment of glaucoma.

PRODUCT NAME	SKU
20G Straight Laser Probe	620200
23G Straight Laser Probe	620201
25G Straight Laser Probe	620202
27G Straight Laser Probe	620203
23G Flexible Curved Laser Probe	620210
25G Flexible Curved Laser Probe	620211
27G Flexible Curved Laser Probe	620212
20G Illuminated Straight Laser Probe	620220
23G Illuminated Straight Laser Probe	620221
25G Illuminated Straight Laser Probe	620222
27G Illuminated Straight Laser Probe	620223
23G Illuminated Flexible Curved Laser Probe	620229
25G Illuminated Flexible Curved Laser Probe	620230
Cycloprobe (Cyclophotocoagulation)	620231



## Straight Laser Probe:

Item Number: 620200, 620201, 620202, 620203

- Sterile, Box of 10
- Disposable
- Available in 20ga, 23ga, 25ga, 27ga



## Flexible Curved Laser Probe:

Item Number: 620210, 620211, 620212

- Sterile, Box of 10
- Disposable
- Available in 23ga, 25ga, 27ga



## Illuminated\* Straight Laser Probe:

Item Number: 620220, 620221, 620222, 620223

- Sterile, Box of 10
- Disposable
- Available in 20ga, 23ga, 25ga, 27ga



\*Illumination source is not included or provided by LIGHTMED.



## Illuminated\* Flexible Curved Laser Probe:

Item Number: 620229, 620230

- Sterile, Box of 10
- Disposable
- Available in 23ga, 25ga



\*Illumination source is not included or provided by LIGHTMED.

## Cycloprobe Cyclophotocoagulation Probe:

Item Number: 620231

- Sterile, Box of 5
- Disposable



FDA CE



©2020 LIGHTMED DCA900-01

1130 Calle Cordillera | San Clemente, CA 92673 | USA  
T: 949-218-9555 | F: 949-218-9556 | [sales@lightmed.com](mailto:sales@lightmed.com)  
[www.lightmed.com](http://www.lightmed.com)

